

Oligonucleotide synthesis kits

To order, contact Tedia at sales@tedia.com

CATALOG NO: OP100/OP400

Tedia has partnered with international oligonucleotide synthesis equipment experts to introduce high purity solvent kits designed for installation and calibration of ÄKTA™ oligopilot™ plus and OligoPilot™ 400 synthesizers.

Industry-leading Purification & Packaging

The products, quantities and packaging in Tedia Oligonucleotide synthesis kits are selected to match the recommended solvents for initial calibration of ÄKTA™ oligopilot™ plus and OligoPilot™ 400 synthesizers. The solvents feature Tedia's industry-leading purification and packaging processes, and are backed by decades of solvent production experience supporting oligonucleotide synthesis around the world. Tedia coordinates with the equipment installation schedule to deliver solvents on-site, on-time, anywhere in the world.

- Tedia Oligonucleotide synthesis kits are developed with synthesis equipment experts with optimal products and quantities for ÄKTA™ oligopilot™ plus and OligoPilot™ 400 installation and initial validation
- Tedia Oligonucleotide synthesis kits are the one-stop shop for solvents packaged for use and consolidated for efficient shipment throughout the world
- Tedia Oligonucleotide synthesis kits are easy - one order, one price includes delivery worldwide
- Just like the oligopilot™, Tedia will scale solvent delivery with your team as your oligonucleotide production grows: from bottles to bulk containers to tankers, Tedia has your solvent supply covered
- Customize your Tedia Oligonucleotide synthesis kit to match your chemistry. Change the blend, the bottle or the volume to meet your needs

KEY HIGHLIGHTS

Optimized
Efficient
Easy
Scaleable
Customizable



ÄKTA™ oligopilot™ plus 100

PRODUCT	DESCRIPTION	QUANTITY	TOTAL LITERS
AB1124-154	Acetonitrile	6 bottles	15
AB3941-154	Activator Solution	1 bottle	2.5
CB3976-154	Capping A Solution	1 bottle	2.5
CB3982-925	Capping B1 Solution*	1 bottle	1
CB3983-925	Capping B2 Solution*	1 bottle	1
DB0490-154	Deblock Solution	1 bottle	2.5
DB9901-243	Diethylamine in Acetonitrile	1 bottle	1
OB3994-154	Oxidation Solution	1 bottle	2.5

OligoPilot™ 400

PRODUCT	DESCRIPTION	QUANTITY	TOTAL LITERS
AB1124-154	Acetonitrile	20 bottles	50
AB3941-154	Activator Solution	1 bottle	2.5
CB3976-154	Capping A Solution	1 bottle	2.5
CB3982-925	Capping B1 Solution*	1 bottle	1
CB3983-925	Capping B2 Solution*	1 bottle	1
DB0490-154	Deblock Solution	8 bottles	20
DB9901-243	Diethylamine in Acetonitrile	1 bottle	2.5
OB3994-154	Oxidation Solution	2 bottles	5

* Transfer and mix Capping B1 Solution to Capping B2 Solution bottle prior to use

*ÄKTA™ and oligopilot™ are trademarks of General Electric Company